Title: Guidelines for Submitting Specimens to NLHRS Form #: NHRL-HRS-LABF027-8

Date of issue: 2023-06-23



NATIONAL LABORATORY FOR HIV REFERENCE SERVICES

GUIDELINES FOR SUBMITTING SPECIMENS TO NLHRS

	Submission Checklist
Ø	NLHRS Requisition completed
V	Provided up to date & secure email for Final Reports
Ø	Met specimen requirements for selected assay
Ø	Shipment is packaged according to TDG Regulation
Ø	Shipper and Courier are TDG certified

NOTE: Changes highlighted in grey

1. REQUISITION NLHRS – HIV/HTLV Requisitions must be used for all submissions.

Available on http://open.cnphi-rcrsp.ca/gts/faces/public/laboratory.xhtml?labId=10080&lang=en

2. SPECIMEN REQUIREMENTS

Assay	Specimen Type	Volume	Storage/Shipping
2.1. HIV-1/2 PCR 2.2. HTLV-I/II PCR/Serology 2.3. HTLV-I Quantification	EDTA anticoagulant whole blood only	5mL (minimum 2mL with the exception of babies)	Ambient temperature (received within 48 hours of collection)
2.4. HIV-1 Viral Load	Plasma	minimum 2mL	-70°C
2.5. HIV-2 Viral Load	Plasma	minimum 2mL	-70°C
2.6. HIV-1/2 Serology 2.7. HTLV-I/II INNO-LIA only	Plasma/Serum	minimum 0.6mL	Ship at your storage temperature

Note: Plasma/serum samples that are hemolyzed or lipemic may produce erroneous results for serology testing.

The National Laboratory for HIV Reference Services is accredited to ISO 15189 by the Standards Council of Canada for the specific scopes of accreditation published on www.scc.ca.

OTHER SPECIMENS - DBS, CSF, TISSUE ETC.

- Contact the NLHRS prior to collecting to ensure sample suitability.
- Note: Specimens categorized as "other" are not subject to accreditation by the Standards Council of Canada for the specific scopes of accreditation published on www.scc.ca
- 3. SHIPPING (https:/tc.canada.ca/en/corporate-services/acts-regulations/list-regulations/transportation-dangerous-goods-regulations)

- Address: Specimen Receiving
Attn: Dr. John Kim

National Microbiology Laboratory 1015 Arlington Street, Room J1139

Winnipeg, MB, R3E 3R2

- Ship by Transportation of Dangerous Goods (TDG) certified courier directly to the NLHRS.
- Notify the NLHRS of shipments with the tracking # by fax, phone or email (below).

4. FINAL REPORTS

- Reports will be emailed and/or faxed only. Please provide an up to date and secure fax/email address.

Contact: Lab email: nlhrs-lnsrv@phac-aspc.gc.ca

Lab fax:(204) 318-2221Philip Lacap:(204) 789-6522John Ho:(204) 789-6518